




# Product Catalog



At HR, our mission is to  
provide quality products  
that **positively impact  
people's lives.**





# Table of Contents

<b>About HR</b>	<b>7</b>
<b>Our Brands</b>	<b>9 - 31</b>
HR® Lubricating Jelly	10
Surgilube® Surgical Lubricant	12
EcoVue® Ultrasound Gel	15
Acclimate® Gel Warmer	17
CoverAlls™ Probe Covers	18
AquaFlate™ Sterile Water	19
TruCath™ Catheters	21
MTG Catheters	26
Vitus Ostomy Care	31
<b>Technical Information</b>	<b>34-35</b>
<b>Symbols Glossary</b>	<b>37-39</b>





5 grams

**surgilube®**

**surgical lubricant  
sterile bacteriostatic**

**DO NOT USE IN EYES AND EARS.**

To open break tip off.  
Surgilube® sterile, surgical lubricant may be used where  
indicated. Do not use Surgilube® for any other purpose.  
Surgilube® is a registered trademark of Surgilube, Inc.

# About HR

HR® is a family-owned company that develops, manufactures, and supplies quality products for a variety of markets. We believe in providing quality products and customized customer service to help people be their healthiest version of themselves. HR® is an industry expert in healthcare consumables and offers a broad range of products from medical lubricants and ultrasound gel to catheters, ostomy care and more.





NET WT. 113g (4oz.)

**Sterile** • Water Soluble • Greaseless • Unscented

**H•R<sup>®</sup>**  Lubricating Jelly

NIH WT 1136 (ADZ)

sterile • Water Soluble • Greaseless • Unscented

**H.R.**  <sup>TM</sup>  
Lubricating Jelly

NET WT. 113g (4oz.)

**To OPEN ASEPTICALLY:** Sanitize tube tip and cap prior to opening. Remove cap and peel off seal from tube orifice. Replace cap before using.

**Does not contain Chlorhexidine Gluconate.**

**INGREDIENTS:** WATER, PR  
METHYLCELLULOSE, CARBO  
PROPYLADIPATE, COUMIN

*Sterile* • Water Soluble • Greaseless • Unsc

**L.D.** 

**Sterile** • **Water Soluble** • **Greaseless** • **Unso**

**H.R.**  <sup>TM</sup>  
Lubricating Jelly

NET WT. 113g (4oz.)

**Sterile • Water Soluble • Greaseless • Unsc**




# Our Brands

# HR® Lubricating Jelly

HR® Lubricating Jelly is a sterile, general purpose medical lubricating jelly, that is widely used in hospitals and physicians' offices. The water-soluble, non-irritating, fragrance-free formulation features optimal viscosity that adheres well and spreads evenly. Suitable for most applications that require a lubricant. Available in a range of formats and sizes to suit practitioners' preferences and procedure requirements.

- Water soluble
- Premium Viscosity
- Glycerin Free
- Bacteriostatic
- Chlorhexidine Free
- Latex & Dye Free



[Visit Website](#)





## HR® Lubricating Jelly Flip-Top Tubes

HR® Lubricating Jelly Flip-Top Tubes are sterile and accommodates a broad range of medical procedures where a multi-use product is appropriate. Available in 2 oz and 4 oz.

MFR#	Description	Packaging
201	HR® Lubricating Jelly 4oz. (113g) Flip-Top Tube	12 ea/bx - 6 bx/cs
201-144	HR® Lubricating Jelly 4oz. (113g) Flip-Top Tube	144 ea/cs
203	HR® Lubricating Jelly 2oz. (56.5g) Flip-Top Tube	12 ea/bx - 8 bx/cs
203-192	HR® Lubricating Jelly 2oz. (56.5g) Flip-Top Tube	192 ea/cs
210-192	HR® Lubricating Jelly <b>2X</b> 2oz. (56.5g) Flip-Top Tube	192 ea/cs



## HR® Lubricating Jelly OneShot®

HR® Lubricating Jelly OneShot® are sterile, single-use packets designed to improve the application, minimize waste, and significantly reduce the risk of cross-contamination.

MFR#	Description	Packaging
207	HR® Sterile Lubricating Jelly 3g OneShot®	144 ea/bx - 6 bx/cs
207-1728	HR® Sterile Lubricating Jelly 3g OneShot®	1728 ea/cs
208	HR® Sterile Lubricating Jelly 3g CarePac® - 30 pack OneShot®	30 ea/bx - 48 bx/cs
209	HR® Lubricating Jelly 5g OneShot®	144 ea/bx - 6 bx/cs
209-1728	HR® Lubricating Jelly 5g OneShot®	1728 ea/cs
205	HR® Lubricating Jelly 1.25oz. (36g) OneShot®	48 ea/bx - 6 bx/cs
205-288	HR® Lubricating Jelly 1.25oz. (36g) OneShot®	288 ea/cs
205ST	HR® Lubricating Jelly 1.25 oz. (36g) OneShot® Safe Wrap™	48 ea/bx - 4 bx/cs
211	HR® Lubricating Jelly <b>2X</b> 1.25 oz. (36g) OneShot®	288 ea/cs

## HR® Lubricating Jelly Pre-Filled Syringe

HR® Lubricating Jelly 10ml Pre-Filled Syringe is a sterile single-use product, the oral tip design with easy-on & easy-off ergonomic cap is ideal for catheter use and medical kits.



MFR#	Description	Packaging
206-400	HR® Lubricating Jelly 10ml Pre-Filled Syringe	400 ea/cs


# Surgilube® Surgical Lubricant

Using Surgilube® Sterile Surgical Lubricant for your patients just makes sense. With a unique, carbomer-free formulation, Surgilube® Sterile Surgical Lubricant increases patient satisfaction while giving YOU the satisfaction of focusing on what you do best: providing quality care.

## Our unique formulation features:

- Carbomer-free formulation reduces interference with pap test results
- Chlorhexidine Gluconate to support infection control
- A proprietary sterilization process to maintain viscosity and crystal-clear clarity

When you are looking for a quality medical lubricant, insist on the preferred choice of medical professionals for over 90 years that is never available as a private label brand.



[Visit Website](#)





## Surgilube® Single-Use Foilpacs®

When the foilpac® was released in 1961, its “single use” concept and design forever changed the mindset of the medical community. Today, more than 50 years later, as health systems continue to look for ways to reduce the risk of cross contamination, Surgilube’s single use foilpac® continues to lead the way.

MFR#	Description	Packaging
281020543	Surgilube® 3g foilpac®	144 ea/bx - 12 bx/ctn
281020545	Surgilube® 5g foilpac®	144 ea/bx - 6 bx/ctn
281020557-288	Surgilube® 31g foilpac®	48 ea/bx - 6 bx/cs



## Surgilube® 5g Snap-Off Metal Tube

This 5g Metal Tube with elongated tip allows for precise application. The snap-off tip also makes for easy single-handed use.

MFR#	Description	Packaging
281020555	Surgilube® 5g Snap-Off Metal Tube	48 ea/bx - 3 bx/ctn



## Surgilube® Screw Cap Metal Tubes

This classic style tube is the ideal option when a metal tube is preferred over laminate. The self-piercing cap allows for convenient opening and closing.

MFR#	Description	Packaging
281020502	Surgilube® 2oz. (56.7g) Screw Cap Metal Tube	12 ea/bx - 12 bx/ctn
281020536	Surgilube® 4.25oz. (120.49g) Screw Cap Metal Tube	12 ea/bx - 6 bx/ctn



## Surgilube® Flip-Top Tubes

Our Flip-Top cap conveniently allows the tube to be free standing. The large easy open, ergonomic design also enables single-handed use. Available in 2oz and 4.25oz.

MFR#	Description	Packaging
281020502	Surgilube® 2oz. (56.7g) Screw Cap Metal Tube	12 ea/bx - 12 bx/ctn
281020536	Surgilube® 4.25oz. (120.49g) Screw Cap Metal Tube	12 ea/bx - 6 bx/ctn






# EcoVue® Ultrasound Gel

With a 99% natural formulation free of dyes, parabens and propylene glycol, EcoVue® Ultrasound Gel is gentler on patients and delivers exceptional glide and image quality with less mess and waste. Available in sterile and non-sterile options with a variety of packaging choices, technicians can choose the right amount for the procedure, dramatically reducing waste.

Trust EcoVue® to deliver the quality and performance that today's ultrasound procedures demand.

- Natural formulation
- Exceptional image quality
- No dyes, fragrances, or parabens
- No propylene glycol



[Visit Website](#)



## EcoVue® Ultrasound Gel Single-Use Packets

Our single-use packets are easy to use and specifically designed to help reduce the possibility of cross-contamination associated with reusing bottles. In addition, our single-use packets are also available with a sterile SafeWrap® for an added protection against cross-contamination.

MFR#	Description	Packaging
280	EcoVue® Ultrasound Gel 20g Packet Sterile SafeWrap™	48 ea/bx – 6 bx/cs
280NW	EcoVue® Ultrasound Gel 20g Packet Sterile	100 ea/bx – 4 bx/cs
281	EcoVue® Ultrasound Gel 20g Packet Non-Sterile	100 ea/bx – 4 bx/cs
283	EcoVue® Ultrasound Gel 32g FlexPac® Non-Sterile	100 ea/bx – 4 bx/cs



## EcoVue® Ultrasound Gel Multi-Use Gel

EcoVue® offers a variety of multi-use sizes for diagnostic and therapeutic medical ultrasound procedures and allows technicians can choose the right amount for the procedure.

MFR#	Description	Packaging
284	EcoVue® Ultrasound Gel 60g Flip-Top Tube Non-Sterile	12 ea/bx – 12 bx/cs
285	EcoVue® Ultrasound Gel 250g FlexPac® Non-Sterile	12 ea/bx – 6 bx/cs
286	EcoVue® Ultrasound Gel 250mL (8.5oz.) Bottle Non-Sterile	12 ea/bx – 6 bx/cs



## EcoVue® HV Ultrasound Gel

EcoVue's high quality natural formulation is now available in a new high viscosity option, delivering the 99% natural formulation of EcoVue® with a higher viscosity for procedures that demand a thicker gel with even more staying power. Available in a bottle as well as a variety of sterile and non-sterile single-use packets, EcoVue® HV is the latest innovation focused on delivering the quality and performance to meet the changing needs of today's ultrasound procedures.

MFR#	Description	Packaging
380	EcoVue® HV Ultrasound Gel 20g Packet Sterile SafeWrap™	48 ea/bx – 6 bx/cs
380NW	EcoVue® HV Ultrasound Gel 20g Packet Sterile	100 ea/bx – 4 bx/cs
381	EcoVue® HV Ultrasound Gel 20g Packet Non-Sterile	100 ea/bx – 4 bx/cs
383	EcoVue® HV Ultrasound Gel 32g FlexPac® Non-Sterile	100 ea/bx – 4 bx/cs
381	EcoVue® HV Ultrasound Gel 20g Packet Non-Sterile	12 ea/bx – 6 bx/cs




# Acclimate® Gel Warmer

The Acclimate® gel warmer is specifically designed to warm EcoVue® Ultrasound Gel and is the only warmer designed with single-use packets and standard size flexible packaging in mind. Our newly redesigned gel warmer features a digital display to continuously monitor temperature. Temperature range can be adjusted in 1 degree increments from 95° F to 105° F.



MFR#	Description	Packaging
289F	Acclimate® Gel Warmer - FlexPac® (Digital Display)	1/bx - 6 bx/cs
289B	Acclimate® Gel Warmer - Bottles (Digital Display)	1/bx - 6 bx/cs
289S	Acclimate® Gel Warmer - Single-Use Packets (Digital Display)	1/bx - 6 bx/cs




[Visit Website](#)

# CoverAlls™ Probe Covers

CoverAlls™ Ultrasound Probe Covers provide a protective barrier to prevent cross-contamination and meet infection control guidelines. Featuring a telescopic-fold with a 3D box end, sonographers can easily apply gel and cover the transducer while maintaining a sterile field. Available in 48" and 96" lengths, CoverAlls™ Ultrasound Probe Covers help protect and extend the life of your equipment. Each kit includes a probe cover, sterile EcoVue® HV Ultrasound Gel, and elastic bands.



MFR#	Description	Packaging
39548	CoverAlls™ General Purpose Probe Cover, 5"x48" (13x122cm), Sterile, Polyurethane, Telescopically-folded, 3D box end, with 20g EcoVue® HV Ultrasound Gel	24 ea/bx
39596	CoverAlls™ Intraoperative Probe Cover, 5"x96" (13x240cm), Sterile, Polyurethane, Telescopically-folded, 3D box end, with 20g EcoVue® HV Ultrasound Gel	24 ea/bx
39648	CoverAlls™ General Purpose Probe Cover, 6"x48" (15x122cm), Sterile, Polyurethane, Telescopically-folded, 3D box end, with 20g EcoVue® HV Ultrasound Gel	24 ea/bx
39696	CoverAlls™ Intraoperative Probe Cover, 6"x96" (15x240cm), Sterile, Polyurethane, Telescopically-folded, 3D box end, with 20g EcoVue® HV Ultrasound Gel	24 ea/bx




[Visit Website](#)

# AquaFlate™ Sterile Water

AquaFlate™ Sterile Water Pre-Filled Syringes contain sterile water and allow for easy balloon inflation of a foley indwelling catheter. AquaFlate™ Pre-Filled Syringes are made from clear plastic with measuring markers to ensure the precise application of fluid. Simply connect the luer slip tip to the catheter and use the plunger-rod to deliver the sterile water to inflate the catheter.



MFR#	Description	Packaging
250-400	AquaFlate™ Sterile Water - 10ml Pre-Filled Syringe	400 ea/cs



[Visit Website](#)






# TruCath™ Catheters

TruCath™ products allow your customers to cost-effectively manage their medical needs without sacrificing quality. Our understanding of CMS reimbursement protocols allows us to structure solutions that can help address the challenges of the urological supply market. Because we are strictly a manufacturer, HR never competes with HME suppliers for patient acquisition.

- Designed for maximum comfort and ease of use
- Feature polished eyelets for smooth insertion and removal
- Sterile & Single-use
- Not made with natural rubber latex or DEHP
- Color-coded funnels for easy French size identification
- Available in Male, Female and Pediatric sizes




[Visit Website](#)



## Intermittent Catheters

TruCath™ Intermittent Catheters are designed for maximum comfort and ease of use. They feature polished eyelets for smooth insertion and removal. They are sterile, single-use and are not made with natural latex rubber or DEHP [Di (2-ethylhexyl) phthalate]. Available in Male, Female and Pediatric sizes.

MFR#	Description	Packaging
SC0816	TruCath™ Intermittent Male Catheter – 8 Fr, 16"	30/box - 300/case
SC1016	TruCath™ Intermittent Male Catheter – 10 Fr, 16"	30/box - 300/case
SC1216	TruCath™ Intermittent Male Catheter – 12 Fr, 16"	30/box - 300/case
SC1416	TruCath™ Intermittent Male Catheter – 14 Fr, 16"	30/box - 300/case
SC1616	TruCath™ Intermittent Male Catheter – 16 Fr, 16"	30/box - 300/case
SC1816	TruCath™ Intermittent Male Catheter – 18 Fr, 16"	30/box - 300/case
SC0806	TruCath™ Intermittent Female Catheter – 8 Fr, 7.5"	30/box - 300/case
SC1006	TruCath™ Intermittent Female Catheter – 10 Fr, 7.5"	30/box - 300/case
SC1206	TruCath™ Intermittent Female Catheter – 12 Fr, 7.5"	30/box - 300/case
SC1406	TruCath™ Intermittent Female Catheter – 14 Fr, 7.5"	30/box - 300/case
SC1406NC	TruCath™ Intermittent Female Catheter – 14 Fr, 7.5" (no connector)	30/box - 300/case
SC1606	TruCath™ Intermittent Female Catheter – 16 Fr, 7.5"	30/box - 300/case
SC0810	TruCath™ Intermittent Pediatric Catheter – 8 Fr, 10"	30/box - 300/case
SC1010	TruCath™ Intermittent Pediatric Catheter – 10 Fr, 10"	30/box - 300/case



## Intermittent Coude Catheters

TruCath™ Coude Intermittent Catheters are designed to maneuver through obstructions in the urethra such as an enlarged prostate or urethral stricture. They feature polished eyelets and a coude (curved) tip for smooth insertion and removal. The catheter also contains a marked funnel for optimal catheter tip positioning during catheterization. They are not made with natural latex or DEHP. Catheters are single-use and sterile.

MFR#	Description	Packaging
CC1016	TruCath™ Intermittent Coude Catheter – 10 Fr, 16"	30/box - 300/case
CC1216	TruCath™ Intermittent Coude Catheter – 12 Fr, 16"	30/box - 300/case
CC1416	TruCath™ Intermittent Coude Catheter – 14 Fr, 16"	30/box - 300/case
CC1616	TruCath™ Intermittent Coude Catheter – 16 Fr, 16"	30/box - 300/case
CC1816	TruCath™ Intermittent Coude Catheter – 18 Fr, 16"	30/box - 300/case





## Hydrophilic Coude Catheters

The TruCath™ Hydrophilic Intermittent Coude Catheter is designed to maneuver through obstructions in the urethra such as an enlarged prostate or urethral stricture, while providing maximum comfort during catheterization. The integrated water pouch when burst, activates the catheter coating to create an ultra-slippery surface. The touch-free sleeve helps to minimize cross-contamination and is used to guide the catheter during insertion and removal. Catheters are single-use and sterile.

MFR#	Description	Packaging
HC1016	TruCath™ Hydrophilic Coude Catheter – 10 Fr, 16"	30/box - 300/case
HC1216	TruCath™ Hydrophilic Coude Catheter – 12 Fr, 16"	30/box - 300/case
HC1416	TruCath™ Hydrophilic Coude Catheter – 14 Fr, 16"	30/box - 300/case
HC1616	TruCath™ Hydrophilic Coude Catheter – 16 Fr, 16"	30/box - 300/case



## Hydrophilic Catheters

The TruCath™ Hydrophilic Intermittent Catheter is designed to provide maximum comfort during catheterization. The integrated water pouch when burst, activates the catheter coating to create an ultra-slippery surface. The touch-free sleeve helps to minimize cross-contamination and is used to guide the catheter during insertion and removal. Available in Male, Female and Pediatric sizes. Catheters are single-use and sterile.

MFR#	Description	Packaging
HS1216	TruCath™ Hydrophilic Male Catheter – 12 Fr, 16"	30/box - 300/case
HS1416	TruCath™ Hydrophilic Male Catheter – 14 Fr, 16"	30/box - 300/case
HS1616	TruCath™ Hydrophilic Male Catheter – 16 Fr, 16"	30/box - 300/case
HS1816	TruCath™ Hydrophilic Male Catheter – 18 Fr, 16"	30/box - 300/case
HS1006	TruCath™ Hydrophilic Female Catheter – 10 Fr, 7.5"	30/box - 300/case
HS1206	TruCath™ Hydrophilic Female Catheter – 12 Fr, 7.5"	30/box - 300/case
HS1406	TruCath™ Hydrophilic Female Catheter – 14 Fr, 7.5"	30/box - 300/case
HS1606	TruCath™ Hydrophilic Female Catheter – 16 Fr, 7.5"	30/box - 300/case
HS0810	TruCath™ Hydrophilic Pediatric Catheter – 8 Fr, 10"	30/box - 300/case
HS1010	TruCath™ Hydrophilic Pediatric Catheter – 10 Fr, 10"	30/box - 300/case



## Closed System Catheters

The TruCath™ Closed System Catheter contains a pre-lubricated catheter and attached collection bag for a no-touch catheterization procedure. The introducer tip shields the catheter from bacteria in the beginning of the urethra where most bacteria are concentrated. Sold as a sterile, single-use catheter/bag unit.

MFR#	Description	Packaging
CSC10	TruCath™ Closed System Catheter – 10 Fr, 16"	100/case
CSC12	TruCath™ Closed System Catheter – 12 Fr, 16"	100/case
CSC14	TruCath™ Closed System Catheter – 14 Fr, 16"	100/case
CSC16	TruCath™ Closed System Catheter – 16 Fr, 16"	100/case



## Closed System Catheter Kits

The TruCath™ Closed System Catheter Kit provides the ultimate in convenience to catheterize anywhere at any time. Each kit contains a pre-lubricated catheter and attached collection bag, along with vinyl powder-free gloves, underpad, PVP and BZK wipes – all the supplies needed for a no-touch catheterization procedure while on the go or in the comfort of home.

MFR#	Description	Packaging
CK10	TruCath™ Closed System Catheter Kit – 10 Fr, 16"	100/case
CK12	TruCath™ Closed System Catheter Kit – 12 Fr, 16"	100/case
CK14	TruCath™ Closed System Catheter Kit – 14 Fr, 16"	100/case
CK16	TruCath™ Closed System Catheter Kit – 16 Fr, 16"	100/case



## Insertion Kits

TruCath™ Catheter Insertion Kits are designed for ultimate convenience in a self-catheterization program. Each kit contains the supplies needed to catheterize while on the go or in the comfort of home. They do not include a catheter and are designed for single-use.

MFR#	Description	Packaging
HRIK001	TruCath™ Intermittent Catheter Insertion kit includes vinyl powder-free gloves, 5g lube jelly packet, BZK wipe, underpad and disposal bag	100/case
HRIK002	TruCath™ Intermittent Catheter Insertion kit includes vinyl powder-free gloves, 5g lube jelly packet, BZK wipe, underpad, drainage bag with connector*, and disposal bag	100/case



## Leg & Drainage Bags

TruCath™ Premium leg drainage bags feature comfortable flocked backing, easy to operate flip drainage valve, sterile fluid pathways, and comes with convenient, comfortable, preattached, fabric, elastic leg straps. Available in 500ml and 1,000ml.

TruCath™ 2,000ml night drainage bags are sterile, diamond-shaped, vented drain bags that feature needleless sample ports and a T-tap drainage value. The double hanger is ideal for hanging on a wheelchair, stretcher, or bed. Measurement markings allow for easy determination of urinary output. Latex-free.

MFR#	Description	Packaging
HRLB500	TruCath™ 500ml Premium Flocked Back Leg Bag with Flip Valve	24/case
HRLB1000	TruCath™ 1,000ml Premium Flocked Back Leg Bag with Flip Valve	24/case
HRNB2000	TruCath™ 2,000ml Night Drainage Bag with Double Hanger and Anti-Reflux Valve	20/case



# MTG Catheters


MTG offers unique male, female, and pediatric catheters, specializing in closed system catheters, and patented catheterization products for paraplegics and those with limited dexterity. Many users find MTG catheters to be life-changing, bringing independence to their bladder care program.



## MTG EZ-Gripper®

The perfect catheter for anyone with limited hand dexterity that can be used without help from a caregiver. Our patented gripping device affixed to the outside of the collection bag assists in easy advancement and our EZ-Advancer® valve holds the catheter stable. A privacy bag is included for discreet disposal.

MFR#	Description	Packaging
52112	16" catheter length with kit components – 12 Fr	50/case
52114	16" catheter length with kit components – 14 Fr	50/case
52116	16" catheter length with kit components – 16 Fr	50/case
50112	16" catheter length without kit components – 12 Fr	50/case
50114	16" catheter length without kit components – 14 Fr	50/case
52812	18" catheter length with kit components – 12 Fr	50/case



[Visit Website](#)



## MTG EZ-Advancer®

Our patented MTG EZ-Advancer® locking valve prevents retraction of the catheter during insertion. This ready-to-use closed system is both easy-to-use and UTI protective. Also available in our Mini-Pak version for users on the go.

MFR#	Description	Packaging
31114	16" catheter length, firm, with kit components (PVP swabs) – 14 Fr	100/case
32112	16" catheter length, firm, with kit components (BZK swabs) – 12 Fr	100/case
32114	16" catheter length, firm, with kit components (BZK swabs) – 14 Fr	100/case
32116	16" catheter length, firm, with kit components (BZK swabs) – 16 Fr	100/case
32212	16" catheter length, soft, with kit components (BZK swabs) – 12 Fr	100/case
32214	16" catheter length, soft, with kit components (BZK swabs) – 14 Fr	100/case
30112	16" catheter length, firm, without kit – 12 Fr	100/case
30114	16" catheter length, firm, without kit – 14 Fr	100/case

## MTG EZ-Advancer® Mini-Pak

MFR#	Description	Packaging
42112	16" catheter length, firm, with kit components (BZK wipe) – 12 Fr	100/case
42114	16" catheter length, firm, with kit components (BZK wipe) – 14 Fr	100/case
42116	16" catheter length, firm, with kit components (BZK wipe) – 16 Fr	100/case



## MTG Eagle® Board

The MTG Eagle® Catheter Board is a unique, innovative product that makes it easier for anyone with poor hand dexterity to self-catheterize, such as quadriplegics (tetraplegics), amputees, people with MS or Parkinson's Disease, and others who have difficulty holding their penis and advancing a urinary catheter.

*Designed to work exclusively with the MTG EZ-Gripper® or MTG EZ-Advancer®.*

MFR#	Description	Packaging
100	MTG Eagle® Board with penile sizing guide, foam pads, carry bag, allen wrench and instructions included	1

## MTG Cath-Lean®

The ideal catheter for women; Protection and Discretion. Closed system and introducer tip help reduce the risk of Urinary Tract Infection recurrence. Mini-Pak fits in purse or pocket.



MFR#	Description	Packaging
42412	6" female catheter length, with kit components – 12 Fr	100/case
42414	6" female catheter length, with kit components – 14 Fr	100/case
40412	6" female catheter length, without kit components – 12 Fr	100/case
40414	6" female catheter length, without kit components – 12 Fr	100/case



## MTG Jiffy Cath®

Closed system urinary intermittent catheters, inclusive of all-important introducer tips, are designed to reduce the risk of recurrent urinary tract infections. For travel, school, or work, reduce the concern for contamination because your hand never touches the catheter tube. Available in both straight and coude as well as Mini-Pak version.

MFR#	Description	Packaging
22112	16" catheter length with kit components – 12 Fr	100/case
22114	16" catheter length with kit components – 14 Fr	100/case
22116	16" catheter length with kit components – 16 Fr	100/case
20112	16" catheter length, without kit components – 12 Fr	100/case
20114	16" catheter length, without kit components – 14 Fr	100/case
20116	16" catheter length, without kit components – 16 Fr	100/case

## MTG Jiffy Cath® Straight Tip Mini-Pak

MFR#	Description	Packaging
12112	16" catheter length, with kit components – 12 Fr	100/case
12114	16" catheter length, with kit components – 14 Fr	100/case
12116	16" catheter length, with kit components – 16 Fr	100/case

## MTG Jiffy Cath® Coude Tip

MFR#	Description	Packaging
22612	16" catheter length, coude tip, with kit components – 12 Fr	100/case
22614	16" catheter length, coude tip, with kit components – 14 Fr	100/case
22616	16" catheter length, coude tip, with kit components – 16 Fr	100/case
22712	16" catheter length, coude tip, with kit components – 12 Fr, SOFT	100/case
22714	16" catheter length, coude tip, with kit components – 14 Fr, SOFT	100/case
22716	16" catheter length, coude tip, with kit components – 16 Fr, SOFT	100/case
20612	16" catheter length, coude tip, without kit components – 12 Fr	100/case
20614	16" catheter length, coude tip, without kit components – 14 Fr	100/case





## MTG Kiddie-Kath®

Closed system catheter with introducer tip helps reduce risk of urinary tract infection recurrence. Kids can learn to catheterize themselves quickly, making restroom independence a reality. For wheelchair users, there isn't a need to transfer to the toilet to catheterize.

MFR#	Description	Packaging
32108	10" catheter length, 700ml collection bag, w/kit (BZK wipe) – 8 Fr	100/case
32110	10" catheter length, 700ml collection bag, w/kit (BZK wipe) – 10 Fr	100/case
32508	16" catheter length, 1500ml collection bag, w/kit (BZK swabs) – 8 Fr	100/case
32510	16" catheter length, 1500ml collection bag, w/kit (BZK swabs) – 10 Fr	100/case
30108	10" catheter length, 700ml collection bag, without kit – 8 Fr	100/case
30110	10" catheter length, 700ml collection bag, without kit – 10 Fr	100/case



## MTG Stoma Stoppers

Prevent stoma leakage and reduce Stoma Stenosis with our Stoma Stoppers, a “must have” for everyday stoma care. The stoma stopper comes with a hydrocolloid dressing you seal over the stoma cap.

MFR#	Description	Packaging
ACE08	8 Fr Stoma Stopper	1
ACE10	10 Fr Stoma Stopper	1
ACE12	12 Fr Stoma Stopper	1
ACE14	14 Fr Stoma Stopper	1
AFD30	Hydrocolloid Dressings	30/box



## Aquaflush® Transanal Irrigation

Aquaflush transanal irrigation makes bowel emptying easy. A solution is using an Aquaflush transanal irrigation system that makes it easy to empty your bowel safely using rectal irrigation, also known as transanal irrigation.

[REF] 2810208.23



**sterile bacteriostatic  
surgical lubricant**  
**DO NOT USE IN EVERETT EARS**  
*This product available in other pharmacy.*  
**NET WT. 3g**

[REF] 281020543

**Surgilube®**

**surgical lubricant  
sterile bacteriostatic  
DO NOT USE IN EYES AND EARS**  
*This product available in tubes at pharmacy*  
**NET WT. 3g**

[REF] 281020543



**REF 281020543**  
**surgical lubricant**  
**sterile bacteriostatic**  
**DO NOT USE IN EYES AND EARS**  
*This product available in tubes at pharmacy*  
**NET WT. 3g**

REF 201020543



REF 281020543  
surgical lubricant  
sterile bacteriostatic  
DO NOT USE IN EYES AND EARS  
this product available in tubes at pharmac  
NET WT. 3g

CONTENTS

**surgilube®**

REF: 281020543  
surgical lubricant  
sterile bacteriostatic  
DO NOT USE IN EYES AND EARS  
This product available in tubes at pharmacies  
NET WT. 3g

**08**


**Surgilube®**

**sterile bacteriostatic**  
**DO NOT USE IN EYES AND EARS**  
**ET WT. 3g**

# Vitus Ostomy Care

Vitus Ostomy offers superior quality ostomy care products at an affordable price designed to help reduce leakage and protect your skin.

Ostomates can “Live Life Freely” and enjoy peace of mind in knowing that our ostomy products allow them to cost-effectively manage their stoma without sacrificing quality.



[Visit Website](#)



## Vitus Ostomy Skin Barrier Rings

Vitus Ostomy Conforming Adhesive Seals are ideal for any type of ostomy. Available in three sizes, our skin barrier rings create a tight seal to prevent leaks, even with high liquid output. They are conformable and stretchable without unnecessary mess and residue. A great alternative to barrier paste.

MFR#	Description	Packaging
520	Vitus Ostomy 2" Skin Barrier Rings - Thin	20 ea/bx - 20 bx/cs
521	Vitus Ostomy 2" Skin Barrier Rings - Standard	20 ea/bx - 20 bx/cs
522	Vitus Ostomy 4" Skin Barrier Rings	10 ea/bx - 20 bx/cs



## Vitus Ostomy Crescent Barrier Extensions

Vitus Ostomy Crescent Barrier Extensions are curved adhesive strips designed to prevent the edges of barrier from lifting and help keep your pouch secure. The thin, hydrocolloid formulation helps conform to uneven body curves. Crescent Barrier extensions can be used in place of tape and act as a frame, providing an extra layer of protection to minimize leaks. They are highly absorbent and can be used under or over pouch wafer.

MFR#	Description	Packaging
513	Vitus Ostomy Crescent Barrier Extensions	30 ea/tray - 50 tray/cs



## Vitus Ostomy Barrier Paste

Vitus Ostomy Barrier Paste is pectin-based and serves as a caulk to fill uneven body contours. It can help prevent leakage and skin irritation around the stoma and provide an optimal seal for a better fit and longer wear time, even with high liquid output. Also available in an alcohol free, no sting formulation.

MFR#	Description	Packaging
500	Vitus Ostomy 2oz. Barrier Paste	24 tubes/box - 2 bx/cs
501	Vitus Ostomy 2oz. No Sting Barrier Paste	24 tubes/box - 2 bx/cs





## Vitus Ostomy Lubricating Deodorant

Vitus Ostomy Lubricating Deodorant helps eliminate odors while lubricating the pouch to make emptying and cleaning easier. Can be used with drainable or closed pouches. Available in three convenient sizes.

MFR#	Description	Packaging
580	Vitus Ostomy 8oz. Lubricating Deodorant Bottle	12 ea/bx - 6 bx/cs
581	Vitus Ostomy 8oz. Lubricating Deodorant Spray Bottle	12 ea/cs
582	Vitus Ostomy .267oz/7.9ml Lubricating Deodorant Packet	30 ea/bx - 48 bx/cs



## Vitus Adhesive Remover Wipes

Vitus Adhesive Remover Wipes gently cleans and removes adhesives from the skin. Each box contains 50, individually wrapped 2" x 2" wipes.

MFR#	Description	Packaging
7050R	Vitus Adhesive Remover Wipes	50 ea/bx - 24 bx/cs



## Vitus Protective Barrier Wipes

Vitus Protective Barrier Wipes provide a barrier film layer on skin under tapes, skin barriers, adhesives, and hydrocolloid dressings to help protect against irritation. Each box contains 50, individually wrapped wipes.

MFR#	Description	Packaging
6050B	Vitus Protective Barrier Wipes	50 ea/bx - 24 bx/cs

# Technical Information

### **HR® Lubricating Jelly**

Instructions for Use

Safety Data Sheet

### **Surgilube® Surgical Lubricant**

Instructions for Use

Safety Data Sheet

### **EcoVue® Ultrasound Gel**

Instructions for Use

Safety Data Sheet

### **Acclimate® Gel Warmer**

Operating Instructions

### **CoverAlls™ Probe Covers**

Instructions for Use

### **AquaFlate™ Sterile Water**

Instructions for Use

### **TruCath™ Intermittent Catheters**

Instructions for Use

### **MTG Catheters**

Instructions for Use

### **Vitus Ostomy Care**

Instructions for Use





32g .600"

**EcoVue**  
Ultrasound

NON-STERILE

If irritation occurs, discontinue use & consult a physician.  
Emergo Europe - Philadelphia 20  
2514 AP The Hague - The Netherlands

0085180900283

LOT 78399  
2024-06-30

Rev 3021-1 PATENT PENDING

Symbols & information available at [EcoVue.com](http://EcoVue.com)  
Hill Pharmaceuticals, Inc. • 2020 Eastern Blvd.  
Suite 201 • York, PA 17402 USA  
877-302-1110

**EcoVue**  
Ultrasound Gel

NON-STERILE

If irritation occurs, discontinue use & consult a physician.  
Emergo Europe - Philadelphia 20  
2514 AP The Hague - The Netherlands

0085180900283

LOT 78399  
2024-06-30

Rev 3021-1 PATENT PENDING

Symbols & information available at [EcoVue.com](http://EcoVue.com)  
Hill Pharmaceuticals, Inc. • 2020 Eastern Blvd.  
Suite 201 • York, PA 17402 USA  
877-302-1110

**EcoVue**  
Ultrasound Gel

NON-STERILE

If irritation occurs, discontinue use & consult a physician.  
Emergo Europe - Philadelphia 20  
2514 AP The Hague - The Netherlands

0085180900283

LOT 78399  
2024-06-30



















Rev 3021-1 PATENT PENDING

Symbols & information available at [EcoVue.com](http://EcoVue.com)  
Hill Pharmaceuticals, Inc. • 2020 Eastern Blvd.  
Suite 201 • York, PA 17402 USA  
877-302-1110














# Symbols Glossary

## Symbols Glossary

United States' Title 21 of the Code of Federal Regulations (CFR) Section 801.109(b)(1)			
Symbol	Reference	Title of Symbol	Description
	21 CFR Section 801.109 (b) (1) and United States Food and Drug Administration Guidance for Industry – Alternative to Certain Prescription Device Labeling Requirements, Issued on January 21, 2000	Prescription Use Only	CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner
ISO 15223-1:2016 – Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements			
Symbol	Reference	Title of Symbol	Description
	5.1.1	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC, and 98/79/EC
	5.1.2	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.
		Date of Manufacture	Indicates the date when the medical device was manufactured.
	5.1.4	Use-by date	Indicates the date after which the medical device is not be used.
	5.1.5	Batch number	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	5.1.6	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified
	5.2.1	Sterile	Indicates a medical device that has been subjected to a sterilization process.
	5.2.3	Sterilized by ethylene oxide treatment	Indicates a medical device that has been sterilized using ethylene oxide.
	5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
	5.2.6	Do Not Re-sterilize	Indicates a medical device that is not to be re-sterilized.
	5.2.7	Non-Sterile	Indicates a medical device that has not been subjected to a sterilization process.
	5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	5.3.2	Keep away from sunlight	Indications a medical device that needs protection from light sources.
	5.3.4	Keep Dry	Indicates a medical device needs to be protected from moisture.
	5.3.7	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.

## Symbols Glossary

	5.4.4	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	5.4.5 and Annex B.2 and Guidance for Industry and US Food and Drug Administration Staff - Recommendations for Labeling Medical Products to Inform Users that the Product or Product Container is not Made with Natural Rubber Latex; Issued on December 2, 2014	Not made with natural rubber latex	Indicates that neither the medical device nor the packaging of the medical device contains the presence of natural rubber latex.
<b>European Medical Device Directive 93/42/EEC / European Medical Device Regulation 2017/745</b>			
	Article 17	CE Conformity Marking and Notified Body Number	Product conforms to the applicable requirements for a medical device as set forth in the regulation and assessed by the certifying notified body.
	Article 17	CE Conformity Marking	Product conforms to the applicable requirements for medical device as set forth in the regulation and assessed by the certifying notified body.
	Annex 1, 23.3(a)	Single Sterile Barrier System	Sterile barrier system
	Annex 1, 23.3(a)	Double Sterile Barrier System	Double entry package with two sterile barrier systems
	Annex 1, 23.2(q)	Medical Device	Designation that the item is classified as a Medical Device
<b>BS EN 15986:2011 – Medical devices – Symbol for use in the labelling of medical devices. Requirements for labelling of medical devices containing phthalates</b>			
	BS EN 15986:2011 EN 15986:2011(E) Annex B	Does not contain DEHP	Indicates a medical device that does not contain the phthalate plasticizers DEHP
	BS EN 15986:2011 Reference no. A.5	Contains or presence of phthalate DEHP	Indicates presence of Bis (2-ethylexyl) phthalate (DEHP)
<b>Other Symbols</b>			
	NA	Made in the USA	All or virtually all significant parts and processing of this product are of U.S. origin.
	NA	Australian Sponsor	Indicates the authorized sponsor in the Australian market.



**Corporate Headquarters**

2600 Eastern Blvd #201  
York, PA 17402

**Manufacturing Facility**

3630 Concord Rd  
York, PA 17402

**Get in Touch**

(877) 302-1110 (toll-free)  
(717) 252-1110

[inquiries@hrpharma.com](mailto:inquiries@hrpharma.com)

[HRpharma.com](http://HRpharma.com)